

HP 3396 Integrator II

BASIC Quick Reference Manual

Formatted by Martin Hepperle, 2021

CONTENTS

CONTENTS	1	FORMAT diskspec, volume, number_files	11
CONVENTIONS	4	GET filespec	11
LOG-ON PROCEDURE	4	INET_CONFIGURATION	11
LEAVING BASIC	4	JOIN filespec	11
HELP MESSAGES	4	LOAD filespec	11
HALTING A RUNNING PROGRAM	4	PACK diskspec	11
ABORTING A COMMAND	4	PURGE filespec	11
OPERATORS	4	RENAME filespec1, filespec2	11
NUMERIC VARIABLES	5	R[UN] [{ filespec line_number }]	11
TYPING TEXT	5	SAVE filespec	11
LABELS	5	SET EXT [-] event_number	11
MULTIPLE STATEMENTS ON A LINE	6	SET MARGIN column_number	11
Strings and Substrings	6	SYSTEM	11
NUMERIC ARRAYS	7	TERMREAD	12
STRING ARRAYS	7	DEBUGGING AND RUN CONTROL COMMANDS	12
AUTOCALL PROGRAMS	8	ENTER BREAK COMMAND	12
DEVICE ADDRESSES	8	ABORT	12
FILE SPECIFIERS AND TYPES	8	BREAK { line_number OFF }	12
COMMAND RANGE	9	DEBUG [{ ON OFF }]	12
EDIT AND LIST COMMANDS	9	LET variable_name = value	12
AUTO_NUMBER	9	RESUME	12
DATE [MM/DD/{YY}YY]	9	SET variable_name = value	12
D[ELETE] [range][,range][, ...]	9	STEP [{ ON OFF }]	12
DQ [range][,range][, ...]	9	TRACE { range OFF }	12
EXIT	9	GENERAL PROGRAM STATEMENTS	12
HELP	10	BEEP	12
L[IST] [range][,range][, ...]	10	BREAK	12
MODIFY [line_number]	10	DATA value[, value][, ...]	12
N[EXT]	10	DEBUG { ON OFF }	12
NOLIST	10	DEF func_name (1 to 3 params) = func_definition	13
P[RINT] [expression_list]	10	DIM array_name or string_name (size) [...]	13
PROMPT { ON OFF }	10	END	13
RENUMBER	10	IMAGE: format_string	13
SCRATCH	10	INPUT variable_name_list	13
SERIALIZE	10	INTEGER { variable array } [...]	13
TIME [HH:MM:SS]	10	[LET] variable_name = expression	13
XADDRESS [device_address]	10	LINE INPUT string_variable	13
XLIST [range][,range][, ...]	10	OPTION BASE { 0 1 }	13
XMARGIN [column_number]	10	PRINT [expression] [{ , ; } expression] [...]	13
FILE AND DEVICE COMMANDS	10	PRINT USING { string line_number label }: expression_list	13
COPY filespec1, filespec2	10	Formatting Numbers	14
CREATE filespec, filesize	10	Formatting Strings	14
DIR[ECTORY] [{ diskspec filespec }]	11	RANDOMIZE	14
		RESTORE [line_number label]	14

READ variable_list	15	OPEN #channel: NAME { filespec	
REM [any text] or ! [any text]	15	device_address }	21
SET run_parameter value	15	PACK diskspec	21
SET EXT [-] event_number	15	PRINTER IS # channel	21
SET FONT n	15	PRINT #channel : expression_list	21
SET MARGIN column_number	15	PRINT #channel USING { string line_number	
SET RUNNUM run_number	15	label } : expression_list	21
STOP	15	PURGE filespec	21
TAB (column_number)	15	READ #channel : variable_list	21
TRACE { line_number OFF }	15	RELEASE REST #channel	21
LOOPING AND CONTROL FLOW STATEMENTS	15	RENAME filespec1, filespec2	21
DO	15	SET #channel: MARGIN column_number	21
DO UNTIL test	15	SET #channel: RECORD numeric_expression	21
DO WHILE test	15	PLOTTING STATEMENTS	22
ELSE	15	PLOT X,Y	22
EXIT DO	15	PEN OFF	22
EXIT FOR	16	PEN ON	22
FOR index= first TO last [STEP size]	16	END PLOT	22
GOSUB { line_number label }	16	MATH AND STRING FUNCTIONS	22
GOTO { line_number label }	16	ABS(x)	22
IF test THEN action [ELSE alternative] END IF	16	ANGLE(x,y)	22
LOOP	16	ATN(x)	22
LOOP UNTIL test	16	BSTR\$(x,n)	22
LOOP WHILE test	16	BVAL(A\$,n)	22
ON index GOSUB { line_number label } [...]	16	CHR\$(x)	22
ON index GOTO { line_number label } [, ...]	16	COS(x)	22
RETURN	17	DATE	22
WAIT [DELAY] numeric_expression	17	DATES	23
WAIT TIME { numeric_expression		DATE2\$	23
string_expression }	17	DEVICES\$(x)	23
EXCEPTION PROTECTION STATEMENTS AND		EPS(x)	23
FUNCTIONS	17	EXP(x)	23
CAUSE EXCEPTION exception_number	17	FLOAT(four_character_string)	23
CONTINUE	17	FLOAT\$(x)	23
END EXCEPTION	17	FP(x)	23
END WHEN	17	INT(x)	23
EXLINE [({ line_number label })]	17	INTRND(x)	23
EXTTEXT\$ (exception_number)	17	IP(x)	23
EXTYPE	19	KEY\$	23
RETRY [{ ALL ({ line_number label }) }]	19	LCASE\$(A\$)	23
USE	19	LEN(A\$)	23
WHEN EXCEPTION IN	19	LOG(x)	23
FILE AND DEVICE STATEMENTS	19	LTRIM\$(A\$)	23
ASK #channel: DATUM string_variable	20	MAX(x,y)	23
ASK #channel: LENGTH numeric_variable	20	MAXNUM	23
ASK #channel: LIFTYPE numeric_variable	20	MIN(x,y)	23
ASK #channel: RECORD numeric_variable	20	MOD(x,y)	23
ASSIGN filespec, key_number	20	NUM(A\$)	23
CHAIN filespec	20	ORD(A\$)	23
CLOSE #channel	20	PI	23
COPY filespec1, filespec2	20	POS(A1\$,A2\$)	23
CREATE filespec, filesize	20	REAL(x)	24
DIRECTORY { diskspec filespec }	20	RND	24
ERASE REST #channel	21	ROUND(x,n)	24
HPIL_IO send, receive	21	RTRIM\$(A\$)	24
INET_IO address, command, string_variable	21	SERIAL_NUMBER\$	24
		SGN(x)	24

SIN(x)	24	IDENTIFIERS\$	27
SQR(x)	24	INJTIMES	27
STR\$(x)	24	ISTDAMT	27
TAN(x)	24	ISTDNUM	27
TIME	24	METHOD_NAMES\$	27
TIMES\$	24	MULT	27
UCASE\$(A\$)	24	NAME\$(calibration_number)	27
UND (numeric variable)	24	NUMCALB	27
USING\$(A\$,x)	24	NUMGRPS	27
VAL(A\$)	24	NUMLEV(calibration_number)	27
BINARY FUNCTIONS	24	NUMPEAKS	27
BINAND (x1,x2)	24	PEAKNUM(calibration_number)	27
BINCMP(x)	24	PROC\$	27
BINEOR(x1,x2)	24	PROCFILE\$	27
BINIOR (x1,x2)	24	REPORT_FILES	27
ROTATE(x,n)	24	REPORT_MEMOS\$	27
SHIFT(x,n)	24	RF(peak_number)	27
CHROMATOGRAPHIC STATEMENTS AND FUNCTIONS	25	RT(peak_number)	27
Statements	25	RUNFILE\$	27
ANALYZE filespec [, I]	25	RUNNUM	28
AREA_PERCENT	25	SAMPAMT	28
BTL_BCD	25	SAMPNAME\$	28
BTL_BIN	25	SAMPNUM	28
GETCALIB filespec	25	SEPCODE\$(peak_number)	28
GETMETH filespec	25	TITLE\$	28
GETSEQ filespec	25	UNCALRF	28
REPORT	25	WIDTH(peaknumber)	28
SAVEMETH filespec	25	SIGNAL DATA ACCESS STATEMENTS AND FUNCTIONS	28
SAVESEQ filespec	25	Statements	28
START RUN_LATER	25	INIT_ACCESS #channel	28
START RUN_LATER END key_number	25	INC_SLICE_NUM (numeric_expression)	28
START RUN_NOW	25	INC_SLICE_TIME (numeric_expression)	28
START RUN_NOW END key number	26	Functions	28
START SEQ_LATER	26	SLICE_AREA	28
START SEQ_LATER END key number	26	SLICE_TIME	28
START SEQ_NOW	26	SLICE_WIDTH	28
START SEQ_NOW END keynumber	26	SLICE NUM	28
Functions	26	SIGNAL_LEVEL	28
AMT (peak_number)	26	RUN PARAMETERS	29
AMT_LBL\$	26	Command or Statement	29
AREA (peak_number)	26	SET run_parameter value	29
CALAMT(calibration_number, level_number)	26	Functions or Parameters	29
CALFIT\$	26	AR_REJ	29
CALNUM(peak_number)	26	ATT2	29
CALRF(calibration_number, level_number)	26	CHT_SP	29
CALRT (calibration_number)	26	PK_WD	29
CALTYPE\$ (peak_number)	26	THRSH	29
EXT\$(event_number)	27	ZERO	29
GROUP_NAMES\$(group_number)	27		
GROUP_SUM(group_number)	27		
HEIGHT(peak_number)	27		

CONVENTIONS

[]	Optional terms.	...	More of the same.
<i>Italics</i>	Supply a value.	{ }	Select one of the items (is separator).

LOG-ON PROCEDURE

1. Type **BA** at the system prompt (*). BASIC is now active on the HP 3396 SERIES II keyboard.
2. Type **BX** at the system prompt. BASIC is now active on an external display terminal. See Section 15 of the User Manual for connection instructions. The HP 3396 SERIES II keyboard is inactive until External BASIC is exited, except that CTRL-BREAK may be used to reset the instrument.

LEAVING BASIC

To exit BASIC (either local or external), type **EXIT** or just **E**. If the message **KEEP PROGRAM IN WORKSPACE [Y/*N]** : appears, respond **Y** (or just **ENTER**) to return to system mode without erasing the workspace program (it will be ready for use the next time BASIC is entered) or **N** to clear the workspace and return to system mode.

HELP MESSAGES

The **HELP** reminder (**TYPE "H" FOR HELP**) appears immediately after logon. Help may be requested whenever the BASIC prompt (**>**) or the break mode prompt (**ENTER BREAK COMMAND**) is present by typing **H** or **HELP**.

The HP 3396 SERIES II displays a list of section headings. The first four sections list names used in the program. Note that names no longer used are removed only when **RUN** is executed. The remaining sections contain listings of functionally related BASIC commands, keywords, functions and parameters. The message **ENTER SECTION NUMBER:** follows the list.

Respond with the number of the desired section. Additional information for that section is displayed. The same result is obtained by typing **H section_number**.

HALTING A RUNNING PROGRAM

In local (HP 3396 SERIES II keyboard) BASIC, the **BREAK** key suspends program execution, preserves all variables, and enters break mode. In external BASIC (using an external display terminal), the **BREAK** key is not active on the terminal; use **CTRL-Y** instead. See the **DEBUGGING AND RUN CONTROL** section.

ABORTING A COMMAND

The **BREAK** key halts command execution as soon as safely possible and returns the BASIC prompt **>**. If a command or BASIC line is typed in, but the **BREAK** key is used before **ENTER**, the input is ignored.

OPERATORS

Execution is from left to right for operators at the same priority level. Use parentheses to override the priority order. Nested parentheses are evaluated from the inside out; six levels of nested parentheses

are permitted. Where alternate syntactical forms are shown here or elsewhere in this manual, it means that all forms are accepted as input, but are converted to the preferred form (given first).

In Decreasing Priority Order

\wedge (or **)	Exponentiation.
NOT , unary -	Logical negation, arithmetic negation.
* , / , DIV , MOD	Multiplication, real division, integer division (no remainder), remainder (modulo).
+ , -	Addition, subtraction.
& (or +)	String concatenation.
Relationals	= , # , < , > , <= , >= <> and >< are acceptable alternate forms of “not equal to.”
AND	Logical conjunction.
OR , XOR	Logical disjunction, exclusive OR.

NUMERIC VARIABLES

Names must start with a letter and may contain up to 31 characters, including letters, numbers, and the underscore character. All characters are significant. Lowercase letters are upshifted. Names cannot be the same as certain keywords and functions. Variables are handled as either 16-bit integers (-32768 to 32767) or 32-bit floating point numbers (-1.7E38 to 1.7E38). Conversions between the two forms are handled automatically.

TYPING TEXT

At turn on, the letter keys print in uppercase (capital letters). Holding down the SHIFT key while pressing a letter key prints in lowercase. Pressing CTRL-C (hold CTRL down and press C) reverses the current sense of the SHIFT key for letters, i.e., if an unshifted letter key printed uppercase before CTRL-C was pressed, it will print lowercase afterward. The SHIFT key then selects uppercase. For keys other than letter keys, the SHIFT key always selects the superscripted symbol on the keycap regardless of whether CTRL-C has been used; the lower symbol is printed if the SHIFT key is not used.

LABELS

A label is defined by its appearance as the first or only item in a line, and it must be followed by a colon (or a @ sign). A label name follows the same set of rules that are followed for simple numeric variable names. Labels may be used as transfer destinations before they are defined.

A label can be used instead of a line number as the destination of a statement such as **GOTO** or **GOSUB**. Statement 350 below causes control to be transferred to the label appearing in line 300.

```
300 LABEL1: PRINT X
310 Statement
320 Statement
...
350 GOTO LABEL1
```

When forming a list of destinations, line number and label parameters may be mixed. For example,

```
100 ON X GOTO 100, HERE, 200, THERE
```

MULTIPLE STATEMENTS ON A LINE

Two or more statements may be used on the same program line, provided they are separated by colons (:), semicolons (;), @ signs, or commas (,). There are some exceptions. Note that the restriction on ELSE only applies when it is part of a multi-line IF THEN ELSE construct. A label may be used in front of all statements.

Alone	First on line	Last on line
DATA	DO	ELSE
DIM	DO	WHILE USE
DEF	DO	UNTIL END IF
INTEGER	FOR	LOOP
IMAGE	Labels	NEXT
OPTION BASE	WHEN EXCEPTION IN	REM

Strings and Substrings

Strings are sets of characters enclosed by single or double quote pairs. String length (number of characters) is limited only by the amount of memory available. String variable names consist of up to 31 characters (as for numeric variables), followed by a dollar sign (\$). Unlike numeric variables and arrays, a simple (non-array) string variable and a string array may not have the same name.

All strings must be dimensioned in DIM statements which specify the maximum number of characters, or physical length, of the string. The actual number of characters is the logical length (logical length of A\$, below, is 6). If an assignment exceeds the physical length, the string is not changed and an exception is generated.

```
1 DIM A$(9), ASTRING$(5)
2 LET A$ = "string"
3 ASTRING$ = 'az"by'
```

Strings can be compared by relational operators on an ASCII numerical basis, and concatenated by the & (or +) operator.

Character positions in strings are numbered from leftmost = 1. Substrings (parts of a string) are specified by one or two numbers or expressions in parentheses. The first number specifies a beginning character.

Two numbers separated by a colon specify beginning and ending characters respectively. If there is no second number, the substring extends to the end of the parent string. Note that the colon must be present (line 10).

Two numbers separated by a semicolon specify the beginning character and the number of characters in the substring respectively. If A\$ is the string "12345", then B\$, C\$, and D\$ are identical to "345".

```
10 B$ = A$(3:)
11 C$ = A$(3;5)
12 D$ = A$(3;3)
```

Assume a substring expression SUB\$(X:Y). If X is greater than Y a null string is returned. See line 103. The examples assume the parent string SUB\$="ABCDEFGH". Asterisks enclose the substrings.

```
100 PRINT "*" & SUB$(1:1) & "*" *      *A*
101 PRINT "*" & SUB$(2:5) & "*" *      *BCDE*
102 PRINT "*" & SUB$(2: ) & "*" *      *BCDEFGH*
103 PRINT "*" & SUB$(4:1) & "*" *      **
```

When a substring appears on the left side of a LET assignment, the parent string characters presently in the positions specified by the substring are replaced by the assigned value. If the first subscript is

greater than the second subscript, as in line 203, the value assigned is inserted starting at the beginning position designated by the first subscript. Existing characters are pushed to the right to create space. If the resulting parent string exceeds its dimensioned length, an exception is generated.

When only one position is specified in the substring, characters in the parent string are completely replaced, beginning at the designated position. Note in line 204, that there are more parent string characters to be replaced than there are characters to fill the positions. The extra characters in the parent string, specifically "G" and "H", are deleted.

The following examples all operate on the parent string "ABCDEFGH".

Statement	Resulting SUB\$
201 SUB\$(2:2) = "Q"	"AQCDEFGH"
202 SUB\$(2:2) = "QR"	"AQRDEFGH"
203 SUB\$(3:2) = "X"	"ABXCDEFGH"
204 SUB\$(3:) = "ACUS"	"ABACUS"
205 SUB\$(2:3) = ""	"ADEFHG"
206 SUB\$(2:3) = "123"	"A123DEFGH"
207 SUB\$(2:0) = "12"	"A12BCDEFGH"

NUMERIC ARRAYS

Array storage is allocated by the **DIM** statement (arrays of numbers which may be integer or floating point) or the **INTEGER** statement (arrays of integers). An array name may be the same as a simple (non-array) variable name in the same program and follows the same rules. The first parameter in a dimensioning statement (**DIM** or **INTEGER**) defines the number of rows in the array. Second and third parameters, if present, define the number of columns and number of "layers" (the third dimension) respectively.

```
100 DIM B1(100), BA_OH(2,15)
101 INTEGER E(5,50,2)
```

STRING ARRAYS

String arrays can be one- or two-dimensional. A string array name follows the same form as for a single string. Line 200 defines a one-dimensional array with 2 elements; each element may contain up to 10 characters. Line 210 fills the array from a (not shown) **DATA** statement. In line 220, the subscript (1) specifies the element number to be printed.

```
200 DIM A1$(2)(10)
210 READ A1$(1),A1$(2)
220 PRINT A1$(1)
```

Line 300 defines a two-dimensional array with 3 rows and 2 columns. Each of the 6 elements may contain 5 characters. Lines 310 through 350 fill the array from a (not shown) **DATA** statement. The subscripts that identify the array element must be enclosed in parentheses.

```
300 DIM A2$(3,2)(5)
310 FOR I=1 TO 3
320   FOR J=1 TO 2
330     READ A2$(I,J)
340   NEXT J
350 NEXT I
```

Substrings of elements are specified by one or two additional subscripts in a second set of parentheses.

```
400 SUB1$ = A1$(2)(30)
401 SUB1$ = A1$(2)(6:6)
402 SUB2$ = A2$(2,1)(2:5)
```

AUTOCALL PROGRAMS

A program assigned to key 0 (zero) is an *Autocall* program and will run automatically at the end of each analytical run. *Autocall* programs cannot accept keyboard input and certain statements are not allowed. They are

INPUT	START RUN_NOW [END]
LINE INPUT	START SEQ_LATER [END]
START RUN _LATER [END]	START SEQ_NOW [END]

DEVICE ADDRESSES

Devices attached to the HP 3396 SERIES II are identified by addresses. These are numbers in the range -2 to 30, some of which are reserved (or pre-defined):

Address	Device Type
-2	Raster mode (internal printer)
-1	RS-232C
0	HP 3396 internal printer, alphanumeric
1-7	Reserved for HP-IB devices
8-30	HP-IL loop devices

FILE SPECIFIERS AND TYPES

A complete file specifier (abbreviated filespec) is **D:NAME.EXT**

D: is the disk specifier (abbreviated diskspec), and consists of a disk drive name followed by a colon (:). The default disk specifier used (for commands only) is **M:**, or may be set by typing in the disk drive name followed by a colon. Legal disk drive names are the capital letters **A** through **Z**, with these exceptions: **H** names the Host, **M** names the HP 3396 SERIES II Internal Memory, and **E**, **G**, **I**, **O** and **Q** may not be used. A disk specifier must have a connected device associated with it.

NAME, the name of the file, may contain up to 8 characters if an extension is used or 10 characters without an extension.

EXT is the extensions and must consist of exactly 3 characters. Some commands and keywords add default extensions if none is explicitly provided. Valid extensions, file types, and file descriptions are shown in the following table. "HP 3396 SERIES II" files are in an LIF format which is unique to this instrument.

Ext	Ext ID	LIF ID	LIF ID	File Type
BAA	5_	1	0001	LIF ASCII BASIC program file
BAS	4_	-10241	D7FF	HP 3396 BASIC program file in memory Image format
DAT	9_	-10241	D7FF	HP 3396 BASIC data file
MET	1_	-10241	D7FF	HP 3396 Method file
SEQ	2_	-10241	D7FF	HP 3396 Sequence file
CAL	3_	-10241	D7FF	HP 3396 Calibration file
RAW	6_	-10241	D7FF	HP 3396 Unbunched signal data file
PRO	7_	-10241	D7FF	HP 3396 Processed peak file from an analytical run
PRA	B_	-10241	D7FF	HP 3396 Processed peak file from an ANALYZE command

RPT	8_	-10241	D7FF	LIF ASCII Report file from an analytical run
RPA	D_	-10241	D7FF	LIF ASCII Report file from an ANALYZE command
BNC	A_	-10241	D7FF	HP 3396 Bunched data from an analytical run
BNA	C_	-10241	D7FF	HP 3396 Bunched data from an ANALYZE command
UA1	H_	-10241	D7FF	LIF ASCII User-defined type
UA2	I_	-10241	D7FF	LIF ASCII User-defined type
UA3	J_	-10241	D7FF	LIF ASCII User-defined type
UD1	E_	-10241	D7FF	HP 3396 User-defined type
UD2	F_	-10241	D7FF	HP 3396 User-defined type
UD3	G_	-10241	D7FF	HP 3396 User-defined type

Files with no extension are in LIF ASCII format. LIF binary files may be read but are not created by the HP 3396 SERIES II. The 2-character Extension ID is used instead of the 3-character file name extension in the internal directory structure.

COMMAND RANGE

The term range, which appears as a parameter for several BASIC commands, is defined as

```
line_number [TO | / | -] line_number]
```

It is also possible to use the keywords `F[IRST]` and `L[AST]` to define a range.

Some examples are:

```
list 55
list 20/50
list FIRST TO 60
list F-L
list 1 TO LAST
```

EDIT AND LIST COMMANDS

Commands may be abbreviated by entering enough of the initial letters to distinguish between commands, with the exceptions of `FORMAT`, `NO LIST` and `SERIALIZE`.

A “current line” pointer is kept internally which provides a default line for many of these commands. The `HELP`, `LIST` and `PRINT` commands found in this section can also be used in break mode (see `DEBUGGING AND RUN CONTROL` section).

AUTO_NUMBER

Supply line numbers automatically as program is entered into workspace. A dialog requests the parameters.

DATE [MM/DD/{YY}YY]

With no parameter, lists the current date and time from the HP 3396 SERIES II internal clock. With the parameter, sets the clock to the specified date. The first two digits of the year are optional until the year 2000.

D[DELETE] [range],[range][, ...]

Remove one or more lines from the workspace program.

Default: delete the current line.

DQ [range],[range][, ...]

Delete Quietly one or more lines from the workspace program. Lines removed are not listed. Default: delete the current line.

EXIT

Terminate BASIC and return to system mode.

HELP

See the introductory material in this manual.

L[IST] [range][,range][, ...]

Print all or part of the workspace program on the internal printer. Default: print entire program.

MODIFY [line_number]

List specified line for editing; use **R**, **D**, or **I** to replace, delete, or insert text. Default: edit the current line; insert characters entered. Position left edge of print head below first character to be modified (or to have text inserted in front of). Terminate modify with ENTER.

N[EXT]

Increment the “current line” pointer and list the new current line.

NOLIST

Prevent workspace program from being listed or revised to provide security. Use before **SAVE**; permits program to be saved and run in BAS format, but not listed or revised. Scans program to check for certain errors, so that checking need not be repeated for subsequent runs. Before using **NOLIST**, it is usually desirable to load a copy of the program into workspace using the **LOAD** command to strip off remarks.

P[RINT] [expression_list]

Evaluate and print the expression(s). If no expression is given, list the current line. expressions can consist of or contain BASIC program variables if the variables have allocated data space. Data space is allocated during and after a run until a change fs made in the program or an **EXIT** from BASIC is executed.

```
PRINT (X+Y)/0.01, A$
```

PROMPT { ON | OFF }

Turns printing of BASIC prompt on or off.

RENUMBER

Renumber all or part of the workspace program. A dialog requests parameters. This also allows blocks of lines to be moved.

SCRATCH

Erase the entire workspace program. Same as DQ F/L (delete quiet first to last).

SERIALIZE

Reserved for future use. See the BASIC Reference Manual.

CAUTION: Use of this command without proper hardware may make the program in the BASIC workspace impossible to run.

TIME [HH:MM:SS]

With no parameter, lists the current time from the HP 3396 SERIES II internal clock. With the parameter, sets the clock to the specified time (24 hour format).

XADDRESS [device_address]

Specify the device for **XLIST** and **XMARGIN**. The built-in printer is address -2. Default device: first printer on HP-IL loop.

XLIST [range][,range][, ...]

Print all or part of the workspace program on the device specified by **XADDRESS**. Default: print entire program.

XMARGIN [column_number]

Set right margin on device specified by **XADDRESS**. Default: no margin (infinite line length allowed).

FILE AND DEVICE COMMANDS

COPY filespec1, filespec2

Copy contents of *filespec1* to *filespec2*. *filespec2* is created by **COPY**; an error occurs if it already exists.

CREATE filespec, filesize

Create a file of the size specified (in bytes) on the specified or defaulted disk. An error occurs if file cannot fit on the disk. If *filesize* is

negative, **CREATE** builds a file as large as possible.

DIR[ECTORY] [{ diskspec | filespec }]

Print the directory for all files on a disk (use disk specifier only) or for a specific file (use file specifier). Without parameter the memory disk MDISK (M:) is printed

FORMAT diskspec, volume, number_files

Prepare a disk to record information and assign it a volume name. Maximum length for volume is 6 characters. *number_files* can be any Integer; it is rounded up to the next increment of 8 when formatting floppy or hard disks and to the next increment of approximately 11 for the memory disk (M:)

```
FORMAT A: ,Vo15,84
```

GET filespec

Erase workspace and load specified program file. Valid extensions are BAA and BAS; if none is provided, first BAS and then BAA are tried as defaults. Only the part of a BAA file that fits in available memory is loaded. If a BAS file does not fit, an error occurs.

INET_CONFIGURATION

Prints a table showing the characteristics of the INET devices on the HP-IL/INET loop.

JOIN filespec

Merge a BAA file with the workspace program. If line number duplication occurs, the line from *filespec* takes precedence; non-conflicting lines Interleave.

LOAD filespec

Erase the workspace, load a BAA file and strip of all **REM** or **!** comments and blank lines.

PACK diskspec

Recover fragmented, unused space on a disk.

PURGE filespec

Delete a file from a disk.

RENAME filespec1, filespec2

Change *filespec1* to *filespec2*. An error occurs if *filespec2* already exists. Both files must be on the same disk.

R[UN] [{ filespec | line_number }]

Fetch and execute program *filespec*, if specified. Valid extensions are BAA and BAS. If no extension is given, first BAS and then BAA are tried as defaults. If no file specified, run the workspace program starting at line number. If no line specified, run starting at the first line.

SAVE filespec

Store the workspace program on the default or specified disk. Valid extensions are BAA and BAS. Default extension: BAS.

```
>save M:test
```

```
>save M:TEST.BAA
>dir

VOLUME NAME: MDISK      DRIVE: M
DATE: JAN  1, 1901  00:33:32

FILE NAME LENGTH  CREATED/VERSION
SIG_BUFF.RAW  2048 01/01/01 00:00:00
TEST .BAS  512 01/01/01 00:29:56
TEST .BAA  256 01/01/01 00:33:29

FILES USED FREE MAX
BYTES 2816 75520 80640
```

SET EXT [-] event_number

Set an external event on an HP 19405A S/ECM module on the HP-IL loop to ON or OFF. *event_number* must evaluate to an integer between 1 and 8, inclusive; a negative value indicates OFF state.

SET MARGIN column_number

Define the right margin for printing. Initial default is 80; range is 1 to 256; values >255 give no margin (infinite line length allowed).

SYSTEM

Print a table showing HP 3396 SERIES II configuration, such as available devices on HP-IL.

TERMREAD

After opening a channel with NAME "-1" (-1 is the RS-232C port), may be used to read from the port buffer into A\$ until an end condition occurs.

DEBUGGING AND RUN CONTROL COMMANDS

These commands are valid responses to the break mode prompt

ENTER BREAK COMMAND

The BREAK, DEBUG, STEP and TRACE commands can also be used outside of break mode. The LIST and PRINT listing commands may also be used in break mode, as can the SET run parameter command. Variable values may be listed using the PRINT command to aid debugging.

ABORT

Terminate program execution and return the BASIC prompt.

BREAK { line_number | OFF }

Set a breakpoint at *line_number*. Program suspends just before executing the specified line and break mode is entered.

Only two breakpoints are permitted; if a third one is specified, the older of the two existing breakpoints is cancelled. BREAK OFF cancels all existing breakpoints.

DEBUG { ON | OFF }

DEBUG ON, or just DEBUG alone, activates BREAK and TRACE statements within the program to help isolate programming errors; DEBUG OFF disables them.

LET variable_name = value

Assign a value to a variable. Attempted multiple assignment (LET A=B=5) is treated as

```
LET A=(B=5), which sets A=1 if B=5 and A=0 if B#5.
```

RESUME

Exit break mode and continue program execution.

SET variable_name = value

Assign a value to a variable.

STEP { ON | OFF }

STEP ON (or STEP alone) sets break points at the beginning of every program line; STEP OFF deactivates them.

TRACE { range | OFF }

For the range specified, list each line number as it is executed. Only two ranges are permitted; if a third range is specified, the older of the two existing ranges is cancelled. TRACE OFF cancels all existing ranges.

GENERAL PROGRAM STATEMENTS

BEEP

Emit a beep.

BREAK

Not valid in Autocall programs.

Place a breakpoint in the program. Program suspends at BREAK statement and enters break mode. See DEBUG statement in this section.

```
190 BREAK
```

```
195 IF 1=100 THEN BREAK
```

DATA value[, value][, ...]

Define data list for the READ statements.

```
200 DATA "VALUES", 1, 2, 3
```

DEBUG { ON | OFF }

Not valid in Autocall programs.

DEBUG ON activates BREAK and TRACE statements; DEBUG OFF disables them. Initial state is OFF.

DEF func_name (1 to 3 params) = func_definition

Define a single line numerical function. The parameters are local to the DEF line and do not conflict with program variables of the same name.

```
260 DEF CUB (X) = X
360 DEF CUBES (X, Y, Z) = X+Y+Z
```

DIM array_name or string_name (size) [...]

Reserve storage for numeric arrays, string arrays, and string variables. Initializes numeric values to a detectable undefined state, and strings to null strings.

```
120 DIM A$(10), B(15,20), ARR$(5,10)
(20)
```

END

Terminate program execution and return to the BASIC prompt. Same as STOP. An implied END is present at the end of any BASIC program.

IMAGE: format_string

Template (printout pattern) for PRINT USING and PRINT USING # statements. In an IMAGE statement, the template is an unquoted format string.

```
100 IMAGE: Date of ### %% 19%%
```

INPUT variable_name_list

Not valid in Autocall programs.

Print a ? prompt; accept input from the keyboard and assign it to the variables listed. Multiple input responses on the same line should be separated by commas.

```
170 INPUT A$, B(3,1), C$, X
```

INTEGER { variable | array } [...]

Declare integer variables or reserve space for integer arrays; initialize values to 0.

```
310 INTEGER A, B(6,6)
```

[LET] variable_name = expression

Assign a value to a variable. Attempted multiple assignment (LET A=B=5) is treated as LET A=(B=5), which sets A=1 if B=5 and A=0 if B#5.

```
300 LET A$ = "STRING"
300 X = 100/Y
```

LINE INPUT string_variable

Not valid in Autocall programs.

Accept an entire line of input, including leading and trailing spaces, from the keyboard and assign it to string variable. No prompt is issued.

```
10 LINE INPUT A$
```

OPTION BASE { 0 | 1 }

Set lower limit for array indices. Default state: 1.

PRINT [expression] [{ , | ; } expression] [...]

Evaluate the expression(s) and send output to the current printer. The , separator spaces values apart in 14 character wide fields, while a ; separator places values as close as possible. The ; separator used at the end of a line suppresses carriage return and line feed. Floating point numbers are printed with a decimal point or in exponential format; they are rounded to 6 digits and trailing zeros in the fractional part of the representation are dropped. Numbers which are in integer form are printed without a decimal point.

```
120 PRINT 5; "A STRING", "ANOTHER"
```

PRINT USING { string | line_number | label } : expression_list

Evaluate the expression(s) and send formatted output to the current printer. Format is controlled by a format string, which may be supplied as a quoted string, as a string variable to which the format string has been assigned, or as a line number or label which contains an IMAGE statement. In the following examples, lines 100, 202, and 301 all produce the same printed result.

```
100 PRINT USING "#.###^^^": X/Y
```

```

200 DIM A$(11)
201 LET A$ = "#.###^ ^ ^ ^"
202 PRINT USING A$: X/Y
300 IMAGE: #.###^ ^ ^ ^
301 PRINT USING 300: X/Y

```

Each character in the format string represents one character position in the output line. Characters are either formatting characters (see the table) which control the appearance of items from the list of items being formatted, or literal characters (any not in the table), which are output unchanged.

Format items are sequences of formatting characters, which must be in the order shown (they may be omitted or in some cases repeated):

Justifiers	>	<
Floating characters	+	- \$
Digit places	*	# %
Period	.	
Number signs	#	
Circumflex accent	^	

A format item ends when a left-to-right scan of the format string encounters an out-of-order formatting character, any non-formatting character, or the end of the format string. If the output line exceeds the margin setting for the device, the line is broken into two lines.

Formatting Numbers

Leading characters are generated according to this table:

Floating character		Character generated	
First	Last	non-negative	negative
-	\$	" \$"	"-\$"
\$	-	"\$ "	"\$-"
-	none	" "	"-"
+	\$	"+\$"	"-\$"
\$	+	"\$+"	"\$-"
+	none	"+"	"-"
\$	none	"\$"	"\$-"
none	none	" "	"-"

A format item may contain only one of each type (sign or \$) of floating character. Digit places are represented by the #, and *

characters. Unused leading digit places are replaced by spaces (the # character), zeros (the % character), or asterisks (the * character). Floating characters float to the right over spaces but not over zeros or asterisks.

An integer format item consists of 1 or more of the same type of digit place character. If a decimal point is added (creating a floating point format item), then all digit places to the right of the decimal point must be the # character. Exponential notation is produced by appending four circumflex accents (^ ^ ^ ^) to an integer or floating point format item. The four accents represent the E, + or - sign, and a two digit exponent.

Exponential format prints one digit to the left of the decimal point (if any), so at least one digit place should be present in this position. Numbers are rounded as necessary to fit the specified format. If a number will not fit the specified format, an attempt is made to print the information by moving the decimal point, changing to exponential notation (if not specified), and/or dropping leading \$ or + characters. If this fails, the format item field is filled with asterisks.

Justifiers may be used to mark the start of a new numeric format item, since only one justifier is allowed per format item. The justifier is replaced by the character to its right. Numbers are always right justified.

Formatting Strings

A string may be output using any format item valid for numbers. This allows the format item to be used for both a column heading and the column data. If the format item begins with a justifier, spaces are added on the right (< justifier) or the left (> justifier) to fill the format item width. If there is no justifier, the string is centered with any leftover space placed on the right.

RANDOMIZE

Generate a new seed for the RND function.

RESTORE [line_number | label]

Reset data pointer to the first item in the DATA statement specified by line number or label. If

no line number or label is specified, reset to the first Item in the first **DATA** statement.

READ variable_list

Read values from the **DATA** list and assign them to the variables.

```
130 READ X, Y, Z, A$
```

REM [any text] or ! [any text]

Program remarks and comments; can be substituted for **REM**. **REM** statements may not be used on the same line following certain key words that do not permit trailing statements; however, **!** without a statement separator may be used in these cases.

SET run_parameter value

Set a *run_parameter* to the specified value. *run_parameter* choices are **AR_REJ**, **ATT2**, **CHT_SP**, **PK_WD**, **THRSH** and **ZERO**. See the Run Parameters section for limits.

```
150 SET CHT_SP 5.0
```

SET EXT [-] event_number

Set an external event on an HP 19405A S/ECM module on the HP-IL loop to **ON** or **OFF**. *event_number* must evaluate to an integer between 1 and 8 inclusive; a negative value indicates **OFF** state.

```
100 SET EXT 1
```

SET FONT n

Where *n* is a numeric expression. If *n* equals 0, the font will be small. If *n* is a non-zero, the font will be large.

SET MARGIN column_number

Define the right margin for printing. Initial default is 80; range is 1 to 255; values > 255 give no margin (infinite line length allowed).

SET RUNNUM run_number

Set the run number.

STOP

Terminate program execution and return to the BASIC prompt. Same as **END**.

TAB (column_number)

Used in **PRINT** statements to move print head to specified column.

```
380 PRINT "LEFT"; TAB (33); "RIGHT"
```

TRACE { line_number | OFF }

Not valid in Autocall programs

TRACE *line_number* lists the statement number as it is executed. **TRACE OFF** cancels the trace. See **DEBUG** statement in this section.

LOOPING AND CONTROL FLOW STATEMENTS

DO

DO marks the beginning of a non-indexed loop and a **LOOP**, **LOOP UNTIL**, or **LOOP WHILE** statement marks its end.

DO UNTIL test

DO UNTIL marks the beginning of a non-indexed loop with a test condition and a **LOOP**, **LOOP UNTIL**, or **LOOP WHILE** statement marks its end. If test is true, passes control to the line following the **LOOP**, **LOOP UNTIL**, or **LOOP WHILE** statement. If test is false, proceeds to the immediately following statement.

DO WHILE test

DO WHILE marks the beginning of a non-indexed loop with a test condition and a **LOOP**, **LOOP UNTIL**, or **LOOP WHILE** statement marks its end. If test is false, continues to the line following the **LOOP**, **LOOP UNTIL**, or **LOOP WHILE** statement. If test is true, proceeds to the immediately following statement.

ELSE

Separates the action and alternative parts of an **IF THEN ELSE** statement.

EXIT DO

When executed, transfers control to the line following the **LOOP**, **LOOP UNTIL**, or **LOOP**

WHILE statement for the present **DO**, **DO UNTIL**, or **DO WHILE** loop.

EXIT FOR

When executed, transfers control to the line following the **NEXT** statement for the present **FOR** loop.

FOR index= first TO last [STEP size]

FOR controls the repeated execution of a group of statements bounded by **FOR** and **NEXT** statements which have the same index (a simple numeric variable). Execution is repeated until any **EXIT FOR** condition is satisfied or the value of Index exceeds last. Then control is transferred to the statement immediately following the **NEXT** statement of the **FOR** loop. If Index is omitted in the **NEXT** statement, **NEXT** is paired with the most recent **FOR** that is still unpaired. Default **STEP** size is + 1. last and size are evaluated dynamically with every iteration. They, along with index, may be changed within the loop.

```
10 FOR I=1 TO N
... Statements
80 IF test THEN EXIT FOR
... Statements
50 NEXT I
```

GOSUB { line_number | label }

Transfer unconditionally to subroutine at line number or label.

GOTO { line_number | label }

Unconditional transfer to specified destination.

IF test THEN action [ELSE alternative] END IF

IF starts a conditional transfer or execution. test is evaluated, and if true, the action is executed. If test is false, action is ignored and alternative, if present, is executed. If test is a numeric expression rather than a relation, test is true if the expression is non-zero. action and alternative may be executable statements or transfer destinations (line numbers or labels). A dangling **THEN** or **ELSE** sets up a multi-line **IF**, and a multi-line **IF** must be ended by **ENDIF**. A nested **IF** cannot be on the same line as the enclosing **IF**.

Both actions and alternatives may consist of multiple statements on the same line, with **:** separators, if the action or alternative is the last item in the line. The entire group of statements is either executed or not executed, depending on the condition of test.

```
100 IF A=B THEN 200
150 IF A=B THEN PRINT B ELSE
160 A=-A: PRINT A
170 ENDIF
```

LOOP

Marks the end of a loop started by a **DO**, **DO UNTIL**, or **DO WHILE** statement. When executed, passes control to the **DO**, **DO UNTIL**, or **DO WHILE** statement which began the loop.

LOOP UNTIL test

Marks the end of a loop started by a **DO**, **DO UNTIL**, or **DO WHILE** statement. When executed, passes control to the **DO**, **DO UNTIL**, or **DO WHILE** statement which started the loop if test is false; otherwise transfers to the statement immediately following the **LOOP UNTIL** statement.

LOOP WHILE test

Marks the end of a loop started by a **DO**, **DO UNTIL**, or **DO WHILE** statement. When executed, passes control to the **DO**, **DO UNTIL**, or **DO WHILE** statement which started the loop if test is true; otherwise transfers to the statement immediately following the **LOOP WHILE** statement.

ON index GOSUB { line_number | label } [...]

Indexed transfer to a subroutine. If index is not between 1 and the number of selections, generate an exception.

```
30 ON Y GOSUB 100,200,START
```

ON index GOTO { line_number | label } [, ...]

Indexed transfer to a program line. If index is not between 1 and the number of selections, go to the next line.

```
20 ON X+1 GOTO 200,300,LOCATION
```


RETURN

Marks the end of a subroutine; causes control transfer to the statement immediately following the `GOSUB` which called the subroutine.

WAIT [DELAY]

numeric_expression

Program operation is suspended until numeric expression seconds have elapsed. numeric expression cannot be greater than 86,399 seconds. The `DELAY` keyword is automatically added if omitted.

```
970 WAIT DELAY AMOUNTOFTIME
```

WAIT TIME {

numeric_expression | string_expression }

Suspends program operation until the time specified. numeric_expression is expressed in seconds since the preceding midnight; If greater than one day (86400 seconds) the value is reduced modulo 86400. string expression is specified in 24-hour clock format ("HH:MM:SS"). If the time specified has already occurred in the present day, operation is suspended until that time the next day.

```
200 WAIT TIME 11500  
210 WAIT TIME "05:30:30"
```

EXCEPTION PROTECTION STATEMENTS AND FUNCTIONS

An exception is a deviation from the normal, or expected, program routine that requires special corrective action. Exception protection programming is used to separate the exception from the regular program routine and carry out appropriate actions to process it. Once the exception has been processed, normal program routine resumes when the exception block is properly exited.

CAUSE EXCEPTION

exception_number

Causes the exception indicated by exception number. Used to test exception protection routines.

CONTINUE

In a `USE` block, clear the exception and transfer to the line following the one which caused the exception.

END EXCEPTION

Used in a `USE` block to clear the exception state; must be followed by a `GOTO` to code outside the `USE` block.

END WHEN

Marks the end of the exception protection block which began at `WHEN EXCEPTION IN`. A `GOTO` on the same line may be used to

transfer control to a line other than the one following `END WHEN`.

EXLINE [({ line_number | label })]

With no parameter, return the line number where the last exception occurred, or 0 if none has occurred. With a line number or label parameter, return 1 if the exception occurred in the specified line or 0 if it did not.

EXTEXT\$ (exception_number)

Return exception message corresponding to exception number. If no message, a null string is returned. To print a list of all exception messages, run the following program:

```
10 FOR I=1000 TO 13000  
20 IF EXTEXT$(I)#" THEN  
30 PRINT I;" : "  
40 PRINT EXTEXT$(I)  
50 PRINT  
60 END IF  
70 NEXT
```

Result:

```
1002: OVERFLOW IN EVALUATING NUMERIC EXPRESSION  
1003: OVERFLOW IN EVALUATING NUMERIC SUPPLIED FUNCTION  
1011: OVERFLOW IN INTEGER ASSIGNMENT  
1106: OVERFLOW IN STRING ASSIGNMENT  
2001: SUBSCRIPT OUT OF BOUNDS  
3002: NEGATIVE NUMBER RAISED TO NONINTEGRAL POWER  
3003: ZERO RAISED TO NEGATIVE POWER
```

3004: LOGARITHM OF ZERO OR NEGATIVE NUMBER
 3005: SQUARE ROOT OF NEGATIVE NUMBER
 3008: ATTEMPT TO EVALUATE ANGLE(0,0)
 3101: UNINITIALIZED VARIABLE ACCESSED (OR INVALID CHROMATOGRAPHIC DATA FUNCTION RESULT)
 3102: INVALID (NULL STRING) PARAMETER
 4001: PARAMETER STRING IS NOT A NUMBER
 4002: ARGUMENT OF "CHR\$" OUT OF RANGE
 4003: ARGUMENT OF "ORD" NOT A VALID CHARACTER OR MNEMONIC
 4005: INDEX IN "TAB" <1 OR >32766
 4006: MARGIN SETTING LESS THAN ONE
 4009: ARGUMENT OF "FLOAT" IS INVALID
 4201: FIRST ARGUMENT OF "BVAL" IS ILLEGAL
 4203: FIRST ARGUMENT OF "BSTR\$" IS ILLEGAL
 4204: SECOND ARGUMENT OF "BVAL" OR "BSTR\$" IS NOT AN EVEN NUMBER FROM 2 TO 72
 4401: ARGUMENT OF "SIN", "COS" OR "TAN" OUT OF RANGE
 4901: OUT OF RANGE SET PARAMETER
 4902: REFERENCE TO NONEXISTENT PEAK
 5000: INSUFFICIENT STORAGE AVAILABLE
 5098: "WHEN EXCEPTION" NESTING TOO DEEP
 5099: PARAMETER STORAGE FULL
 7001: CHANNEL NUMBER NOT IN RANGE 0 TO 99
 7002: INVALID OPERATION FOR CHANNEL ZERO
 7003: NONZERO CHANNEL IN "OPEN" ALREADY ACTIVE
 7004: INACTIVE CHANNEL IN FILE STATEMENT OTHER THAN "OPEN"
 7009: ALL 3 I/O CHANNELS IN USE
 7101: FILE NOT FOUND
 7102: INSUFFICIENT DISC SPACE
 7103: DUPLICATE FILE NAME
 7104: DIRECTORY FULL
 7105: END OF FILE (ON READ)
 7106: DISC DIRECTORY BAD (LENGTH =0)
 7107: FOREIGN OR UNFORMATTED DISC
 7108: READ ERROR (USUALLY CRC)
 7109: NO DISC IN DRIVE
 7110: DISC WRITE PROTECTED
 7111: DISC TIMEOUT
 7112: ATTEMPTED WRITE PAST END OF FILE
 7113: BUFFER SIZE TOO SMALL (FOR PACK)
 7114: NEW DISC INSERTED
 7115: DISC NOT USABLE (TOO MANY BAD SECTORS)
 7116: DISC NOT FORMATTED
 7117: DIRECTORY SIZE TOO LARGE (FOR FORMAT)
 7118: SECTOR OUTSIDE FILE BOUNDS
 7119: INCORRECT NUMBER OF BYTES READ
 7120: DEVICE DOES NOT TALK
 7121: DISC ERROR
 7130: BAD LOOP FRAME
 7131: BAD OUTPUT FRAME
 7132: NO INPUT BUFFER
 7133: NO OUTPUT
 7134: "ETE" ERROR
 7135: I/O DEVICE DOWN
 7136: CANNOT USE/FIND I/O DEVICE
 7137: HPIB I/O DEVICE DOWN
 7140: NO HPIL DEVICES
 7141: HPIL DOWN
 7142: TOO MANY HPIL OR HPIB DEVICES
 7143: HPIL CONFIGURATION ERROR
 7144: DISC CONFIGURATION CHANGED
 7150: INVALID PARAMETER(S)
 7151: CAN NOT GET FILE CONTROL BLOCK
 7152: FILE ALREADY OPEN
 7153: DISC BEING ACCESSED
 7154: ILLEGAL RENAME OPERATION (DIFFERENT DISCS)
 7155: ILLEGAL VOLUME NAME
 7156: NO DISC NAME SPECIFIED
 7157: ILLEGAL DISC NAME
 7158: DISC DOES NOT EXIST
 7159: NO FILE NAME SPECIFIED
 7160: ILLEGAL FILE NAME
 7161: ILLEGAL EXTENSION NAME
 7162: WEAR INDICATOR WARNING
 7163: HOST PROTOCOL ERROR
 7164: HOST CONVERSATION BUSY
 7165: HOST CONVERSATION ABORTED
 7166: ILLEGAL HOST OPERATION
 7167: INVALID DEVICE ADDRESS
 7168: EPROM POWER FAIL
 7170: DEVICE UNAVAILABLE
 7171: HOST UNAVAILABLE
 7180: HOST FILE NOT OPEN
 8001: "READ" BEYOND END OF DATA
 8101: INVALID DATUM FOR "READ" OF NUMBER
 8109: INVALID DATUM FOR "READ" OF STRING
 8201: INVALID FORMAT-STRING
 8202: NO FORMAT-ITEM IN FORMAT-STRING FOR OUTPUT LIST

```

8210: NUMERIC FORMAT-ITEM IS OVER 25 CHARACTERS
9001: POWER FAIL OR LOOP BREAK
9003: LOOP IS DOWN
9005: ADDRESS IS NOT ON LOOP OR DOES NOT SUPPORT KO DATA PATH
9006: ADDRESS IS INACTIVE ON KO DATA PATH
9010: RESULT STRING IS LARGER THAN SPACE ALLOWED
9020: obviously raised when PEEK access a forbidden memory address
9021: INVALID EVENT NUMBER
9025: LOOP IS DOWN
9027: NO EXTERNAL EVENTS INSTRUMENT IN LOOP
9028: EXTERNAL EVENTS INSTRUMENT NOT ACTIVE
9201: INVALID FILE OR I/O DEVICE NAME
9202: BAD .BAS FILE FORMAT
9203: INCOMPATIBLE VERSIONS OF .BAS FILE AND BASIC
9301: HOST HANDSHAKE TIMEOUT
9302: HOST CONVERSATION STAGE FAULT
9303: EXCESS CHARACTERS FROM HOST (BUFFER OVERFLOW)
9305: ILLEGAL CHARACTER RECEIVED FROM HOST
9307: EXCESS AND ILLEGAL CHARACTERS RECEIVED FROM HOST
9309: INVALID MUTE HOST CONVERSATION
9803: END OF SIGNAL
9804: END OF SIGNAL (ABORTED SOURCE)
9805: INCORRECT FILE CONTENT
9806: REQUESTED NONPOSITIVE SLICE NUMBER
9807: SLICE ACCESS NOT INITIALIZED
9808: SLICE AREA OUT OF RANGE
9900: ERROR IN GETTING CAL, METH , OR SEQ FILE
9902: INVALID IN AUTOCALL
9903: INVALID FROM BA/BX
9904: KEY 0 ASSIGNED
9911: RDY CHK ABORTED
9912: EQUIB DELAY ABORTED
9913: SO PATH OFF
9914: METH GET FAILED
9915: LOOP ERROR
9916: STOP PRESS DURING SEQ
10001: INDEX OUT OF RANGE IN ON-GOSUB
10002: RETURN WITHOUT CORRESPONDING GOSUB
10100: RETRY WITHOUT EXCEPTION
10101: "USE" OR "END WHEN" WITHOUT EXCEPTION
12004: ILLEGAL NUMERIC VALUE SPECIFIED FOR TIME-EXPRESSION
12005: ILLEGAL STRING VALUE SPECIFIED FOR TIME-EXPRESSION

```

EXTYPE

Return the exception number. The function `EXTYPE` returns 0 if used outside of an exception processing block.

RETRY [{ ALL | ({ line_number | label }) }]

In a `USE` block, clears the exception and transfers control: with no parameter, to the line which caused the exception; with the `ALL` parameter, to the `WHEN EXCEPTION IN` statement to retry the entire block; with a line number or label parameter, to the specified line, which must be in the protected block.

USE

Separates the block of statements being protected against exceptions from the

correction statements. If no exception has occurred when `USE` is reached, control transfers to the line following the `END WHEN` statement.

WHEN EXCEPTION IN

Starts a block of statements for which exception protection is desired. When an exception occurs, statements following `USE` are executed to perform corrective action.

```

10 WHEN EXCEPTION IN
... Statement
50 USE
... Exception processing statements
90 END WHEN

```

FILE AND DEVICE STATEMENTS

channel may range from 0 to 99. Channel 0 is defined to be the device from which BASIC is being run; channel 0 is always open.

File specifiers, which appear as parameters in some of these statements, may be supplied as either quoted strings or string variables. A complete *filespec* is required; the *diskspec* portion cannot be defaulted.

LIF ASCII files (BAA, RPT, RPA, UA1, UA2, UA3, and files without extensions) and LIF binary files are read a record at a time; if the record is assigned to a numeric variable, conversion to a number is done if possible; otherwise an exception is generated. If the record read is assigned to a string variable but will not fit, an exception is also generated.

The size of elements in DAT files is 4 bytes for all numbers; strings use 4 bytes plus the number of characters in the string.

Checks are done to see that requests to read numbers or strings from DAT files are reading an item of the correct type.

Other HP 3396 SERIES II unique files are read *n* bytes at a time by reading a string or substring of length *n*. No checks are done regarding the format of what is read. It is also possible to read a 4 byte (floating point internal format) number directly from such a file by Reading into a numeric variable.

ASK #channel: DATUM string_variable

Assign the type of the next data item to be read from file *channel* to string variable. For DAT files, returns "STRING", "NUMERIC", "NONE" (end of file reached), or "UNKNOWN" (bad file format or incorrect record pointer position in file); for all other files and devices, returns "UNKNOWN".

```
40 ASK #5: DATUM A$
```

ASK #channel: LENGTH numeric_variable

Assign the physical length in bytes of file *channel* to the variable. If not a file, return 0.

ASK #channel: LIFTYPE numeric_variable

Assign the LIF file type number in file *channel* to the given variable (LIF ASCII: 1; LIF binary: -2; LIF HP 3396 SERIES II unique: -10241_d (D7FF_h); if not a file, 0 is returned).

ASK #channel: RECORD numeric_variable

Assign the current value of the record number pointer in file *channel* to the variable. This is the byte position in the file rather than a true record number. If not a file, return 0. For LIF ASCII (or LIF binary) files, note that the value returned includes bytes used for record size counts and padding, and so will not be the same as the total number of characters in the records that have been read.

```
40 ASK #6: RECORD R
```

ASSIGN filespec, key_number

Assign the program *filespec* to the key *key_number*. *filespec* is a string variable or quoted string; *key_number* ranges from 0 to 9. The program is loaded and run when the key is pressed while in system mode. The program assigned to key 0 is loaded and run after each analytical run.

CHAIN filespec

Terminate execution of currently running program and clear the workspace; load and execute a program file. Only BAS files may be chained.

CLOSE #channel

Close access to file *channel*. If channel is zero or not in use, no action is taken.

COPY filespec1, filespec2

Copies the first file, creating the second file.

CREATE filespec, filesize

Create a file of *filesize* bytes on the specified or default disk. An error occurs if the file will not fit. If *filesize* is negative, CREATE builds the largest possible file.

DIRECTORY { diskspec | filespec }

Print a directory of all files on a disk (use *diskspec* only) or the entry for a specific file (use complete *filespec*).

ERASE REST #channel

Set the end-of-file marker to the current position of the file pointer, effectively erasing the remainder of the file.

HPIL_IO send, receive

Reserved [reverse engineered by MH]

Send one or more HP-IL frames to the loop and receive the result (if any). *Send* must be a string composed of byte pairs encoding the frame payload in the first byte and the type of the frame in the second byte. The frame types are DAB=0 or 1, END=2 or 3, CMD=4, RDY=5. The *receive* parameter must be a string variable and receives the result from the loop. Its first byte is a status byte, which is usually 0x04.

INET_IO address, command, string_variable

Transmit *command* (string or string variable) to device *address* via INET. Response is returned in *string_variable*.

```
20 INET_IO 8, "BOTTLE 100", A$
```

OPEN #channel: NAME { filespec | device_address }

Assign a channel number to a file or an external device. If a file, open it.

```
20 OPEN #20: NAME D$
21 OPEN #30: NAME "B:AFILE.BAS"
22 OPEN #40: NAME "8"
```

PACK diskspec

Recover fragmented, unused space on a disk. *diskspec* must be a string.

PRINTER IS # channel

Define current output device for **PRINT** statements.

PRINT #channel : expression_list

Write data to channel specified.

```
530 PRINT #12: B$, C$, D
```

PRINT #channel USING { string | line_number | label } : expression_list

Evaluate the expression(s) and send formatted output to the channel. If expression list

contains more than one item and channel is a file, all formatted outputs are combined into a single record. Format is controlled by a format string, which may be supplied as a quoted string, as a string variable to which the format string has been assigned, or as a line number or label which contains an **IMAGE** statement. See the **PRINT USING** statement for the format string details.

```
120 PRINT #5 USING "%":A
```

PURGE filespec

Delete a file from a disk.

READ #channel : variable_list

Read data from the specified *channel*. Assign values to the variables listed.

```
120 READ #24: X, Y, Z
```

RELEASE REST #channel

Return unused space at the end of the file (past the end-of-file marker) to the disk. For LIF ASCII files, the file pointer must be at the end-of-file (EOF) marker, or no action is taken. Even if the file is empty, a read must be done to return an EOF error before the system is aware of the EOF position and can execute this command; no action is taken if these conditions are not met.

RENAME filespec1, filespec2

Change *filespec1* to *filespec2*. An error occurs if *filespec2* already exists.

SET #channel: MARGIN column_number

Set right margin for channel. See **SET MARGIN** for range. If this statement is not used, the default is no margin (infinite line length allowed).

SET #channel: RECORD numeric_expression

Set record pointer in file channel to record number *numeric_expression*. Record number is a byte position in the file, rather than a true record number. Use **RECORD 0** to set to the first byte. For LIF ASCII (or LIF binary) files, any non-zero value will set the file pointer to the physical end of the file; it is not possible to

set the file pointer to a particular record except by reading all previous records.

PLOTTING STATEMENTS

The chromatogram plotting operation which is invoked by the **PLOT** key may also be used from BASIC to plot any desired data. This operation starts when a **PLOT X,Y** statement is executed and terminates when **END PLOT** is executed. While plotting, data, all carriage returns and line feeds are stripped from printed material and compressed width printing is used; the print head returns to the curve after each string.

PLOT X,Y

If not already plotting, begin plotting operation and turn **PEN ON**. Advance paper by X units; move the plotter print head to Y. The X range is 0 to 16500 units (1 unit = 1/660 cm¹). Y range is 0 to 1312 units². There are certain physical limits on the rapidity with which X and Y may change and yield accurate plots; see the BASIC Reference Manual for more information and specific limitations.

¹ Which is 66 dots per mm or 1676 dpi, thus filling a paper height of 250 mm

² For a paper width of 210 mm this yields about 6 dots per mm or 152 dpi.

PEN OFF

Turn the plot dot **OFF** on the HP 3396 SERIES II plotter.

PEN ON

Turn the plot dot **ON** on the HP 3396 SERIES II plotter.

END PLOT

End plotting operation.

MATH AND STRING FUNCTIONS

Each function returns a value. **x** and **y** represent numeric expressions; **n** represents a numeric expression which is rounded to the nearest integer (if necessary) before use. All returned strings are shown between double quotes to mark the start and end; these quote characters are not part of the returned string.

ABS(x)

Absolute value of **x**.

ANGLE(x,y)

Angle (in radians) between positive x-axis and a vector joining the origin to the point (**x,y**).

ANGLE(5,5) returns 0.785398

ATN(x)

Arc tangent (in radians) of **x**, i.e., value of the angle whose tangent is **x**.

BSTR\$(x,n)

String representation of **x** using base **n**. (**n** must evaluate to an even integer from 2 to 72.)

BSTR\$(3,2) returns "11"

BVAL(A\$,n)

Value of A\$ interpreted using base **n**. (**n** must evaluate to an even Integer from 2 to 72.)

BVAL("1F",16) returns 31.

CHR\$(x)

ASCII character equivalent of **x**. See Table A.

CHR\$(90) returns "Z".

COS(x)

Cosine of x radians.

DATE

Current Julian date from HP 3396 SERIES II internal clock, in YYDDD format.

DATE returns 85001 (for January 1, 1985).

DATE\$

Current date, in YYYYMMDD format.

DATE\$ returns "19610911" (for September 11, 1961).

DATE2\$

Current date, in MMM DD, YYYY format.

DATE2\$ returns "SEP 11, 1961".

DEVICE\$(x)

Parameter x is an HP-IL address from 8 to 30; a string is returned showing the characteristics of the device. Format is the same as in SYSTEM command listings. If no device is at that address, returns the null string ("").

EPS(x)

Resolution error for x .

EXP(x)

The number e raised to the x power.

FLOAT(four_character_string)

Convert a four_character_string representing a floating point number in Internal format (as may be read from a processed peak file) Into a BASIC number. Reverse of FLOAT(x).

FLOAT\$(x)

Convert a floating point number into a number in four character string representing in a floating point number Internal format. Reverse of FLOAT(four_character_string).

FP(x)

Fractional part of x .

INT(x)

Largest Integer less than or equal to x .

INTRND(x)

Nearest Integer to x .

IP(x)

Integer part (by truncating at the decimal point) of x .

KEY\$

Returns one character for the last key hit on the HP 3396 SERIES II keyboard (or the terminal keyboard in External BASIC mode), then

clears to zero (the NUL character). See Table A.

LCASE\$(A\$)

A\$ with each uppercase letter shifted to lowercase.

LCASE\$("LOWERCASE") returns "lowercase".

LEN(A\$)

Number of characters (logical length) in A\$.

LOG(x)

Natural (base e) logarithm of x .

LTRIM\$(A\$)

A\$ with leading spaces deleted.

MAX(x,y)

Larger (more positive) of x and y .

MAXNUM

Largest floating-point value that can be represented (1. 70141 E+38).

MIN(x,y)

Smaller (more negative) of x and y .

MOD(x,y)

x modulo y ; $x - y * INT(x/y)$

NUM(A\$)

ASCII decimal equivalent of first character in A\$. See Table A.

NUM("&#") returns 42; NUM("&#") returns 42.

ORD(A\$)

Ordinal position of A\$. A\$ must be only 1 character, or else consist of one of the ORD mnemonics in Table A. For example, the ordinal position for backspace is 8 and the mnemonic is "BS".

ORD ("BS") returns a 8.

PI

Represented as 3.14159.

POS(A1\$,A2\$)

Character position of the first character of the first occurrence in A1\$ of the string value A2\$.

If the value of `A2$` is not present in `A1$`, a zero is returned.

`POS("SUBSTRING", "STRING")` returns 3.

REAL(x)

Converts an integer number `x` into the real (floating point) format.

RND

Next number in a pseudo-random sequence.

ROUND(x,n)

`x` rounded to `n` places to the right of the decimal point, or `-n` places to the left if `n < 0`.

RTRIM\$(A\$)

`A$` with trailing spaces deleted.

SERIAL_NUMBER\$

Returns a 1 a-character serial number if an EPROM with a serial number is installed; otherwise returns a null string.

SGN(x)

Algebraic sign of `x`. Negative = -1; 0 = 0; positive = +1.

SIN(x)

Sine of `x` radians.

SQR(x)

Positive square-root of `x`.

BINARY FUNCTIONS

These operate on the Internal 16-bit 2's complement integer representation. MSB is bit 15; LSB is bit 0; there is no carry bit. Floating point numbers. are first converted to the nearest integer. For `SHIFT` and `ROTATE` operations, positive `+n` are shifting/rotating to the right.

BINAND (x1,x2)

Bit-by-bit logical-and of `x1` and `x2`.

BINCMP(x)

Bit-by-bit complement of `x`.

BINEOR(x1,x2)

Bit-by-bit exclusive-or of `x1` and `x2`.

BINIOR (x1,x2)

Bit-by-bit inclusive-or of `x1` and `x2`.

STR\$(x)

String equivalent to `x` as it would be printed from a simple PRINT statement.

`STR$(10)` returns "10".

TAN(x)

Tangent of `x` radians.

TIME

Current time, in seconds since the previous midnight. Resolution is 0.05 seconds.

TIME\$

Current time, in 24-hour clock format (HH:MM:SS).

UCASE\$(A\$)

`A$` with each lowercase letter in uppercase.

`UCASE$("upper")` returns "UPPER"

UND (numeric variable)

1 if variable value is undefined, otherwise 0 is returned.

USING\$(A\$,x)

Formatted representation of number `x`, using `A$` as the template. See `PRINT USING` statement.

VAL(A\$)

Number represented by `A$`. `VAL("20")` returns 20.

ROTATE(x,n)

Rotate 16-bit binary representation of `x`, `n` bit positions; performed with wraparound.

`ROTATE(9, -2)` returns 36

SHIFT(x,n)

Shift 16-bit binary representation of `x`, `n` bit positions; performed without wraparound.

`SHIFT(9, 2)` returns 2

CHROMATOGRAPHIC STATEMENTS AND FUNCTIONS

Statements

ANALYZE filespec [, I]

Not valid in Autocall programs.

Reintegrates raw or bunched run data in *filespec*, which must be a complete file specifier. If *I* is used, uses the peak width profile from the original data storage for reintegration. *filespec* must be in a single string variable or quoted string.

AREA_PERCENT

Generate Area% or Height% report (depending on mode selected in Option 4 of Report Options), using data and parameters from the current method.

BTL_BCD

Returns the bottle number {0-99} as read in the Binary Coded Decimal format from the SAMPLE connector on the rear of the integrator.

BTL_BIN

Returns the bottle number (0-255) as read in the binary format from the SAMPLE connector on the rear of the integrator.

GETCALIB filespec

Erase the current calibration and load the specified calibration file. “.ext” type is automatically added if not present. *filespec* must be a quoted string.

GETMETH filespec

Erase the current method and load the specified method file. “.ext” type is automatically added if not present. *filespec* must be a quoted string.

GETSEQ filespec

Erase the current sequence and load the specified sequence file. “.ext” type is automatically added if not present. *filespec* must be a quoted string.

REPORT

Generate a report using the calculation procedure specified by the currently active method.

SAVEMETH filespec

Store the current method on the device and file specified. *filespec* must be a string variable or quoted string.

SAVESEQ filespec

Store the current sequence on the device and file specified. *filespec* must be a string variable or quoted string.

START RUN_LATER

Not valid in Autocall programs

Programmatic initiation of an analytical run. Run begins when [START] is pressed, and statements following **START RUN_LATER** are suspended until the run ends.

START RUN_LATER END key_number

Not valid in Autocall programs

Same as **START RUN_LATER** except: The program containing this statement must be assigned to a numeric key and run by pressing the numeric key in system mode. When the statement is executed, the program containing it is cleared from the workspace. When the run ends, the program assigned to key number is loaded and executed. The *key_number* program is cleared from the workspace when it ends.

START RUN_NOW

Not valid in Autocall programs.

Programmatic initiation of an analytical run. Run starts immediately, and statements following **START RUN NOW** are suspended until the run ends.

START RUN_NOW END key number

Not valid in Autocall programs.

Same as **START RUN_NOW** except: The program containing this statement must be assigned to a numeric key and run by pressing the numeric key in system mode. When the statement is executed, the program containing it is cleared from the workspace. When the run ends, the program assigned to key number is loaded and executed. The key number program is cleared from the workspace when it ends.

START SEQ_LATER

Not valid in Autocall programs.

Programmatic initiation of a sequence of analytical runs. Sequence begins when [START] is pressed and statements following **START SEQ_LATER** are suspended until the sequence ends.

START SEQ_LATER END key number

Not valid in Autocall programs.

Same as **START SEQ_LATER** except: The program containing this statement must be assigned to a numeric key and run by pressing the numeric key in system mode. When the

Functions

AMT (peak_number)

Concentration or amount determined by the calculation procedure.

AMT_LBL\$

Amount label from Option 4 (Report Options) dialog. Maximum length: 10 characters.

AREA (peak_number)

Baseline-corrected area, in 1/8 microvolt-seconds, of the specified peak.

CALAMT(calibration_number, level_number)

Amount of the calibration entry.

CALFIT\$

Type of calibration fit. Returns "P" (point-to-point), "L" (linear regression), or "N" (non-linear i.e., quadratic).

statement is executed the program containing it is cleared from the workspace. When the sequence ends, the program assigned to key number is loaded and executed. The key number program is cleared from the workspace when the sequence ends.

START SEQ_NOW

Not valid in Autocall programs.

Programmatic initiation of a sequence of analytical runs. Sequence starts immediately, and program statements following **START SEQ_NOW** are suspended until the sequence ends.

START SEQ_NOW END keynumber

Not valid in Autocall programs.

Same as **START SEQ_NOW** except: The program containing this statement must be assigned to a numeric key and run by pressing the numeric key in system mode. When the statement is executed the program containing it is cleared from the workspace. When the sequence ends, the program assigned to key number is loaded and executed. The key number program is cleared from the workspace when the sequence ends.

CALNUM(peak_number)

Calibration table number of the specified processed peak. (Peaks are numbered in order of retention time.) This is the inverse of the **PEAKNUM** function.

CALRF(calibration_number, level_number)

Response factor at specified calibration point.

CALRT (calibration_number)

Calibration table retention time of the calibrated peak.

CALTYPE\$ (peak_number)

Type of calibration peak. Returns "R" (Reference Peak), "S" (Internal Standard Peak), "&" (Reference and Internal Standard Peak), or Blank character (other).

EXT\$(event_number)

Indicates state (“ON” or “OFF”) of the event specified. *event_number* must round to an integer from 1 to 8.

GROUP_NAME\$(group_number)

Name of the specified group in the calibration table. Maximum length: 16 characters.

GROUP_SUM(group_number)

Total of the calculated amounts of the peaks in the specified group.

HEIGHT(peak_number)

Baseline-corrected height of the specified peak, in 1/8 microvolts.

IDENTIFIER\$

Instrument identifier. Set using the system mode **IDENTIFIER** command. Maximum length: 12 characters.

INJTIME\$

Date and time when sample was injected. Maximum length: 22 characters.

ISTDAMT

Amount of the internal standard peak.

ISTDNUM

Calibration number of the internal standard peak.

METHOD_NAME\$

Most recently reported method file name (file name of the active method used in the most recent run). Maximum length: 14 characters. Returns null string if the active method is not stored in a file anywhere.

MULT

Value of the constant factor by which all calculation procedure results are multiplied (**MULTIPLIER** parameter).

NAME\$(calibration_number)

Name of the calibrated peak. Maximum length: 16 characters.

NUMCALB

Maximum calibration number used in the present calibration table.

NUMGRPS

Total number of peak groupings in a sample.

NUMLEV(calibration_number)

Number of levels in the calibration for the calibrated peak.

NUMPEAKS

Number of peaks in the processed data file.

PEAKNUM(calibration_number)

The peak number corresponding to the specified calibration number. This is the inverse of the **CALNUM** function.

PROC\$

Calculation procedure in the current calibration table. Possible returns are **NORM**, **ESTD**, **ESTD%**, **ISTD**, and **ISTD%**, followed by either **-AREA** or **-HEIGHT** to show what measurement was used. If current method is not calibrated, returns **AREA%**.

PROCFILE\$

Name of the processed data file. Maximum length: 14 characters.

REPORT_FILE\$

Name of the report file. Maximum length: 14 characters.

REPORT_MEMO\$

Maximum length: 126 characters. Can come from option 7 default report memo or from sample information table.

RF(peak_number)

Response factor of the specified peak.

RT(peak_number)

Measured retention time of the specified peak, in minutes.

RUNFILE\$

Name of the unbunched or bunched data file. Maximum length: 14 characters.

RUNNUM

Current run number.

SAMPAMT

Amount of the sample (used with external and internal standard calculations) .

SAMPNAME\$

Name of the sample. Maximum length: 12 characters.

SAMPNUM

Number of bottle most recently injected or BCD number from the HP 3396 SERIES II back panel or from the signal data file after reintegration.

SEPCODE\$(peak_number)

Peak Separation Code for specified peak. Maximum length: 4 characters.

TITLE\$

Report title from Report Options dialog. Maximum length: 42 characters.

UNCALRF

Response factor used for uncalibrated peaks.

WIDTH(peaknumber)

Width of the specified peak, in minutes.

SIGNAL DATA ACCESS STATEMENTS AND FUNCTIONS

Statements

INIT_ACCESS #channel

Set the stored slice access pointer to 1. Use after OPENing the file but before accessing stored data.

INC_SLICE_NUM (numeric_expression)

Add numeric expression to the stored slice access pointer. Cause an exception if the pointed-to slice does not exist.

INC_SLICE_TIME (numeric_expression)

Add (numeric_expression/20) seconds to the current SLICE_TIME. If the new time does not fall within or on the end of an existing slice, cause an exception. If the new time is valid, adjust the stored slice access pointer to point to the slice which contains the new time.

Functions

SLICE_AREA

Area, in microvolt-seconds, of the slice pointed to by the stored slice access pointer.

SLICE_TIME

Time after start of run, in 1/20 seconds, at the end of the slice pointed to by the stored slice access pointer.

SLICE_WIDTH

Width, in 1/20 seconds, of the slice pointed to by the stored slice access pointer.

SLICE_NUM

Current value of the stored slice access pointer.

SIGNAL_LEVEL

Voltage level of real-time input signal to HP 3396 SERIES II in millivolts. See the BASIC Reference Manual for more information and specific limitations.

RUN PARAMETERS

The Run Parameter functions return the current values of the run parameter settings. The current values can be changed using the **SET** command (which may also be used as a program statement) in this section.

Command or Statement

SET run_parameter value

Set a run_parameter to the specified value. run_parameter choices are **AR_REJ**, **ATT2**, **CHT_SP**, **PK_WD**, **THRSH** and **ZERO**.

Functions or Parameters

AR_REJ

Area rejection level. Range is from 0 to 2147483647 area counts.

ATT2

Chart attenuation in powers of 2, from -8 to 36.

CHT_SP

Chart speed from 0 to 30.0 cm/min for source or filtered plots, and from 0 to 3.0 mm/peak for Unigram plot.

PK_WD

Peak width setting, from 0.01 to 2.5 minutes.

THRSH

Threshold in powers of 2, from -6 to 28.

ZERO

Chart zero offset from -6 for the left margin to 100 for the right margin.

ASCII Decimal	Graphic Display	ORD Mnemonic	HP3396 Key
0		NUL	CTRL + @
1		SOH	CTRL + A
2		STX	CTRL + B
3		ETX	CTRL + C
4		EOT	CTRL + D
5		ENQ	CTRL + E
6		ACK	CTRL + F
7		BEL	CTRL + G
8	(pad)	BS	CTRL + H
9		HT	CTRL + I
10		LF	CTRL + J
11		VT	CTRL + K
12		FF	CTRL + L
13		CR	ENTER
14		SO	CTRL + N
15		SI	CTRL + O
16		DLE	CTRL + P
17		DC1	CTRL + Q
18		DC2	CTRL + R
19		DC3	CTRL + S
20		DC4	CTRL + T
21		NAK	CTRL + U
22		SYN	CTRL + V
23		ETB	CTRL + W
24		CAN	CTRL + X

ASCII Decimal	Graphic Display	ORD Mnemonic	HP3396 Key
25	BREAK	EM	CTRL + Y
26		SUB	CTRL + Z
27	<=	ESC	ESC
28		FS	
29		GS	
30		RS	CTRL + ^
31		MS	CTRL + _
32	(space)	SP	SPACE
<hr/>			
225			LIST
226			DEL
227			ZERO
228			ATT2
229			CHT SP
230			AR REJ
231			THRSH
232			PK WD
233			EXT()
234			INTG()
235			STOP
236			START
237			PLOT
238			TIME
<hr/>			
240			PREP

ASCII Decimal	Graphic Display	ORD Mnemonic	HP3396 Key
241			EDIT
242			STORE
243			LOAD
244			SEQ
245			METH
246			CALIB

ASCII Decimal	Graphic Display	ORD Mnemonic	HP3396 Key
247			REPORT
248			AREA%
249			OP()

TABLE A. HP 3396 Key and Character Set